



Serving Healthcare Since 1971



# GENERAL SURGICAL INSTRUMENTS CATALOG

[www.volmed.co.uk](http://www.volmed.co.uk)

## About Us

With over 50 years of industry experience, VOLMED is a family-owned company specializing in the supply of high-quality surgical instruments to the medical and surgical markets across the United Kingdom, Europe, and globally. Our extensive expertise enables us to assist with project research and development, ensuring that we meet the evolving needs of our clients.

Our philosophy is simple: we treat our customers as we would like to be treated, which places customer satisfaction at the forefront of our priorities. This approach naturally aligns with our commitment to delivering products of the highest quality. Ultimately, our goal is to empower you to provide the best possible care and achieve improved health outcomes by streamlining the process of stocking your practice.

We take pride in the breadth of our product offerings, understanding that sourcing medical tools and surgical instruments from multiple suppliers can be both time-consuming and challenging. To address this, we strive to be your one-stop shop for all your requirements.

## Our History

Our story began over 50 years ago within a family deeply rooted in the medical instrument manufacturing industry. The foundation of our thriving business was built by a team of highly ambitious professionals with expertise in producing plastic surgery, dental and surgical instruments.

Inspired by the potential for excellence, our founder established VOLMED, a company dedicated to providing outstanding customer service and high-quality products. Since then, our team of dedicated professionals has worked tirelessly to ensure that customers worldwide receive only the finest surgical instruments and related services. We never compromise on quality, continually investing in our team, resources, equipment, and processes. Every detail is meticulously examined to ensure we can adapt to changing demands while maintaining the highest standards of quality.

## Our Promise

At VOLMED, our commitment to excellence is unwavering. We strive to achieve the highest standards in everything we do, working relentlessly to provide exceptional service to our customers and everyone associated with our company. We view it as our responsibility to contribute to a better future for the entire community, and we are dedicated to making this vision a reality.

## Quality

Our commitment to quality is reflected in our CE and UKCA-marked instruments, manufactured in accordance with stringent global quality assurance standards, including ISO 9001:2015, ISO 13485:2016, and the Medical Device Directive 93/42/EEC.

## Our Mission

We are committed to manufacturing and supplying products that meet the highest standards of strength, accuracy, and reliability. Recognizing the critical role our equipment plays in saving lives and improving health, we take this responsibility seriously. Our mission is to simplify the procurement process by offering a comprehensive range of high-quality surgical products under one roof.

## Products

At VOLMED, surgical instruments are at the core of our offerings. We provide a wide range of high-quality surgical products suitable for use in operating rooms, wards, veterinary practices, dentistry, and chiropody. Our instruments are crafted to strict standards using a blend of modern technology and traditional craftsmanship, in full compliance with the Medical Devices Directive 93/42/EEC. Complete traceability is maintained throughout the manufacturing process.

To ensure exceptional product performance and reliability, all our products undergo rigorous engineering analysis and stringent inspection methods. VOLMED instruments are guaranteed for five years against manufacturing defects and are made from premium surgical-grade stainless steel, in accordance with ISO 7153-1:2016 specifications. Explore our product listings and catalogues to discover our comprehensive range of surgical instruments.



Esophagoscopy Tracheotomy

**VOLMED**

**VOLMED**

Serving Healthcare Since 1971


Esophagoscopy  
Tracheotomy

Serving Healthcare Since 1971



34  
26


31,26






Shaft length

|             |           |           |           |           |           |           |           |
|-------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| 30 cm / 12" | V-034-440 | V-034-470 | V-034-500 | V-034-530 | V-034-570 | V-034-610 | V-034-650 |
| 40 cm / 16" | V-034-450 | V-034-480 | V-034-510 | V-034-540 | V-034-580 | V-034-620 | V-034-660 |
| 50 cm / 20" | V-034-460 | V-034-490 | V-034-520 | V-034-550 | V-034-590 | V-034-630 | V-034-670 |
| 60 cm / 24" | V-034-461 | V-034-491 | V-034-521 | V-034-560 | V-034-600 | V-034-640 | V-034-680 |



Chevalier-Jackson




Huber  
V-034-690  
Universal handle



Bruenings

|           |             |
|-----------|-------------|
| V-034-700 | 20 cm / 8"  |
| V-034-710 | 30 cm / 12" |
| V-034-720 | 40 cm / 16" |
| V-034-730 | 50 cm / 20" |
| V-034-740 | 60 cm / 24" |

Interchangeable tubes



Bruenings  
V-034-760  
20 cm / 8"  
Interchangeable tube

# Esophagoscopy


Exchangeable Tips

Auswechselbare Zangenaufsätze


Pinzas intercambiables




Fraenkel  
V-034-770



Baracani  
V-034-780




Seiffert  
V-034-790



Fraenkel  
V-034-800



Scheinmann  
V-034-810



Scheinmann  
V-034-820



V-034-830

V-034-840

Bruenings



V-034-850

V-034-860

Bruenings




Lange  
V-034-870



V-034-880

V-034-890

Krause



Krause  
V-034-900

V-034-910




V-034-920

V-034-930

V-034-940

Cordes



Schumacher  
V-034-950  
With drawing cut



V-034-960

V-034-970


Cordes



Cordes  
V-034-980




Rosenberg  
V-034-990



|            |                   |
|------------|-------------------|
| V-034-1011 | Fig. 00 = 8 mm Ø  |
| V-034-1012 | Fig. 0 = 10 mm Ø  |
| V-034-1013 | Fig. 1 = 12 mm Ø  |
| V-034-1014 | Fig. 2 = 14 mm Ø  |
| V-034-1015 | Fig. 3 = 16 mm Ø  |
| V-034-1016 | Fig. 4 = 18 mm Ø  |
| V-034-1017 | Fig. 5 = 20 mm Ø  |
| V-034-1018 | Fig. 6 = 22 mm Ø  |
| V-034-1019 | Fig. 7 = 24 mm Ø  |
| V-034-1020 | Fig. 8 = 26 mm Ø  |
| V-034-1021 | Fig. 9 = 28 mm Ø  |
| V-034-1022 | Fig. 10 = 30 mm Ø |


Mirror without handle

13.5 cm / 5 1/4"  
Plane, without handle




|            |                   |
|------------|-------------------|
| V-034-1026 | Fig. 00 = 8 mm Ø  |
| V-034-1027 | Fig. 0 = 10 mm Ø  |
| V-034-1028 | Fig. 1 = 12 mm Ø  |
| V-034-1029 | Fig. 2 = 14 mm Ø  |
| V-034-1030 | Fig. 3 = 16 mm Ø  |
| V-034-1031 | Fig. 4 = 18 mm Ø  |
| V-034-1032 | Fig. 5 = 20 mm Ø  |
| V-034-1033 | Fig. 6 = 22 mm Ø  |
| V-034-1034 | Fig. 7 = 24 mm Ø  |
| V-034-1035 | Fig. 8 = 26 mm Ø  |
| V-034-1036 | Fig. 9 = 28 mm Ø  |
| V-034-1037 | Fig. 10 = 30 mm Ø |


Mirror with handle




18.5 cm / 7 1/4"  
Plane, with fixed handle  
Rhodium handle




Michel  
25 cm / 10"  
Moveable mirror




V-034-1039  
12 mm Ø  
Fig. 1



V-034-1041  
14 mm Ø  
Fig. 2




V-034-1042  
16 mm Ø  
Fig. 3



V-034-1043  
18 mm Ø  
Fig. 4



V-034-1024  
10 cm / 4"




V-034-1025  
10 cm / 4"

# Esophagoscopy


Laryngeal Polypus Forceps, Epiglottis Retractor

Kehlkopfpolypenzangen, Deckelhalter


Pinzas para polipos laringeos, depresor para epiglottis




**Kahler**  
**V-034-1260**  
40 cm / 16"  
Tip 4 mm




**V-034-1270**  
30 cm / 12"  
Tip 7 mm




**Papilloma**  
**V-034-1265**  
30 cm / 12"  
Tip 7 mm



**Jackson**  
**V-034-1275** 30 cm / 12", 2 mm  
Straight  
**V-034-1276** 30 cm / 12", 4 mm  
Curved




**Jackson**  
**V-034-1280**  
30 cm / 12", 4 mm  
Angled




**V-034-1244** 12 cm / 4 1/4", 2 mm  
**V-034-1245** 14 cm / 5 1/2", 2 mm  
**V-034-1246** 16 cm / 6 1/4", 2 mm  
**V-034-1247** 18 cm / 7", 3 mm  
**V-034-1248** 20 cm / 8", 3 mm  
**V-034-1249** 30 cm / 12", 4 mm  
**V-034-1250** 40 cm / 16", 4 mm  
**V-034-1251** 50 cm / 20", 4 mm




**Jackson**  
24 cm / 9 1/2"





**V-034-1044**  
1 x 2 Teeth




**V-034-1045**  
Without teeth






**Fraenkel**  
**V-034-1046**  
23 cm / 9"




**V-034-1047**




**V-034-1048**



**V-034-1049**



**V-034-1051**




**Jurasz**  
23 cm / 9"

# Tracheotomy

Trachea Trocars and Dilators, Trachea Retractors

Tracheal-Trokare und Dilatatoren, Tracheal Hakchen


Traqueo tomos Y Dilatadores, Separadores prala traquea




V-026-001  
8 mm

V-026-010  
10 mm

Uckermann  
13 cm / 5"  
Hollow handle




Trousseau  
V-026-020  
14 cm / 5½"




Laborde  
V-026-030  
14 cm / 5½"



Bose  
V-026-040  
7 cm / 2¾"



V-026-050  
18 cm / 7"




V-026-060  
18 cm / 7"




Iterson  
16 cm / 6¼"




V-026-070  
Sharp




V-026-080  
Blunt



Rose  
V-026-090  
13 cm / 5"



Luer  
V-026-100  
11 cm / 4 1/4"




Jackson  
V-026-110  
15 cm / 6"

# Tracheotomy

Tracheal Retractors and Hooks

Tracheal Hakchen und Haken


Seperadores prala traquea Y Ganchos traqueales



Jackson

V-026-120

18 cm / 7"




Lukens

V-026-130

17 cm / 6 3/4"

Blunt




Jackson

V-026-140

17 cm / 6 3/4"

### International scale


|                      |                                                                                   |                                                                                   |                                                                                   |                                                                                   |                                                                                   |                                                                                   |                                                                                    |                                                                                     |                                                                                     |                                                                                     |
|----------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|                      |  |  |  |  |  |  |  |  |  |  |
| Diameter             | 4 mm Ø                                                                            | 5 mm Ø                                                                            | 6 mm Ø                                                                            | 7 mm Ø                                                                            | 8 mm Ø                                                                            | 9 mm Ø                                                                            | 10 mm Ø                                                                            | 11 mm Ø                                                                             | 12 mm Ø                                                                             | 13 mm Ø                                                                             |
| Length of outer tube | 45 mm                                                                             | 50 mm                                                                             | 55 mm                                                                             | 60 mm                                                                             | 65 mm                                                                             | 70 mm                                                                             | 75 mm                                                                              | 80 mm                                                                               | 85 mm                                                                               | 90 mm                                                                               |
|                      | Fig. 0                                                                            | Fig. 1                                                                            | Fig. 2                                                                            | Fig. 3                                                                            | Fig. 4                                                                            | Fig. 5                                                                            | Fig. 6                                                                             | Fig. 7                                                                              | Fig. 8                                                                              | Fig. 9                                                                              |



### Chevalier-Jackson

Outer tube closed  
Inner tube closed

|           | Fig. | Ø     |
|-----------|------|-------|
| V-026-150 | 0    | 4 mm  |
| V-026-160 | 1    | 5 mm  |
| V-026-170 | 2    | 6 mm  |
| V-026-180 | 3    | 7 mm  |
| V-026-190 | 4    | 8 mm  |
| V-026-200 | 5    | 9 mm  |
| V-026-210 | 6    | 10 mm |
| V-026-220 | 7    | 11 mm |
| V-026-230 | 8    | 12 mm |
| V-026-231 | 9    | 13 mm |




### Luer

Outer tube sieve  
Inner tube sieve

|           | Fig. | Ø     |
|-----------|------|-------|
| V-026-240 | 0    | 4 mm  |
| V-026-250 | 1    | 5 mm  |
| V-026-260 | 2    | 6 mm  |
| V-026-270 | 3    | 7 mm  |
| V-026-280 | 4    | 8 mm  |
| V-026-290 | 5    | 9 mm  |
| V-026-300 | 6    | 10 mm |
| V-026-310 | 7    | 11 mm |
| V-026-320 | 8    | 12 mm |
| V-026-321 | 9    | 13 mm |

International scale

|                      |                                                                                   |                                                                                   |                                                                                   |                                                                                   |                                                                                   |                                                                                   |                                                                                   |                                                                                     |                                                                                     |                                                                                     |
|----------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|                      |  |  |  |  |  |  |  |  |  |  |
| Diameter             | 4 mm Ø                                                                            | 5 mm Ø                                                                            | 6 mm Ø                                                                            | 7 mm Ø                                                                            | 8 mm Ø                                                                            | 9 mm Ø                                                                            | 10 mm Ø                                                                           | 11 mm Ø                                                                             | 12 mm Ø                                                                             | 13 mm Ø                                                                             |
| Length of outer tube | 45 mm<br>Fig. 0                                                                   | 50 mm<br>Fig. 1                                                                   | 55 mm<br>Fig. 2                                                                   | 60 mm<br>Fig. 3                                                                   | 65 mm<br>Fig. 4                                                                   | 70 mm<br>Fig. 5                                                                   | 75 mm<br>Fig. 6                                                                   | 80 mm<br>Fig. 7                                                                     | 85 mm<br>Fig. 8                                                                     | 90 mm<br>Fig. 9                                                                     |




Krishaber

|           | Fig. | Ø     |
|-----------|------|-------|
| V-026-330 | 0    | 4 mm  |
| V-026-340 | 1    | 5 mm  |
| V-026-350 | 2    | 6 mm  |
| V-026-360 | 3    | 7 mm  |
| V-026-370 | 4    | 8 mm  |
| V-026-380 | 5    | 9 mm  |
| V-026-390 | 6    | 10 mm |
| V-026-400 | 7    | 11 mm |
| V-026-410 | 8    | 12 mm |



Krishaber-Luer

|           | Fig. | Ø     |
|-----------|------|-------|
| V-026-420 | 0    | 4 mm  |
| V-026-430 | 1    | 5 mm  |
| V-026-440 | 2    | 6 mm  |
| V-026-450 | 3    | 7 mm  |
| V-026-460 | 4    | 8 mm  |
| V-026-470 | 5    | 9 mm  |
| V-026-480 | 6    | 10 mm |
| V-026-490 | 7    | 11 mm |
| V-026-500 | 8    | 12 mm |



Parker

|           | Fig. | Ø     |
|-----------|------|-------|
| V-026-510 | 1    | 16 mm |
| V-026-520 | 2    | 18 mm |
| V-026-530 | 3    | 20 mm |
| V-026-540 | 4    | 22 mm |
| V-026-550 | 5    | 24 mm |
| V-026-560 | 6    | 26 mm |
| V-026-570 | 7    | 28 mm |
| V-026-580 | 8    | 30 mm |
| V-026-590 | 9    | 32 mm |

### International scale

|                      |                                                                                   |                                                                                   |                                                                                   |                                                                                   |                                                                                   |                                                                                   |                                                                                   |                                                                                     |                                                                                     |                                                                                     |
|----------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|                      |  |  |  |  |  |  |  |  |  |  |
| Diameter             | 4 mm Ø                                                                            | 5 mm Ø                                                                            | 6 mm Ø                                                                            | 7 mm Ø                                                                            | 8 mm Ø                                                                            | 9 mm Ø                                                                            | 10 mm Ø                                                                           | 11 mm Ø                                                                             | 12 mm Ø                                                                             | 13 mm Ø                                                                             |
| Length of outer tube | 45 mm                                                                             | 50 mm                                                                             | 55 mm                                                                             | 60 mm                                                                             | 65 mm                                                                             | 70 mm                                                                             | 75 mm                                                                             | 80 mm                                                                               | 85 mm                                                                               | 90 mm                                                                               |
|                      | Fig. 0                                                                            | Fig. 1                                                                            | Fig. 2                                                                            | Fig. 3                                                                            | Fig. 4                                                                            | Fig. 5                                                                            | Fig. 6                                                                            | Fig. 7                                                                              | Fig. 8                                                                              | Fig. 9                                                                              |




Collin

|           | Fig. | Ø     |
|-----------|------|-------|
| V-026-600 | 0    | 4 mm  |
| V-026-610 | 1    | 5 mm  |
| V-026-620 | 2    | 6 mm  |
| V-026-630 | 3    | 7 mm  |
| V-026-640 | 4    | 8 mm  |
| V-026-650 | 5    | 9 mm  |
| V-026-660 | 6    | 10 mm |
| V-026-670 | 7    | 11 mm |
| V-026-680 | 8    | 12 mm |



Fuller

|           | Fig. | Ø     |
|-----------|------|-------|
| V-026-690 | 1    | 18 mm |
| V-026-700 | 2    | 20 mm |
| V-026-710 | 3    | 22 mm |
| V-026-720 | 4    | 24 mm |
| V-026-730 | 5    | 26 mm |
| V-026-740 | 6    | 28 mm |
| V-026-750 | 7    | 30 mm |
| V-026-760 | 8    | 32 mm |



Cubley

|           | Fig. | Ø     |
|-----------|------|-------|
| V-026-770 | 1    | 16 mm |
| V-026-780 | 2    | 18 mm |
| V-026-790 | 3    | 20 mm |
| V-026-800 | 4    | 22 mm |
| V-026-810 | 5    | 24 mm |
| V-026-820 | 6    | 26 mm |
| V-026-830 | 7    | 28 mm |
| V-026-840 | 8    | 30 mm |
| V-026-850 | 9    | 32 mm |



Serving Healthcare Since 1971

## CONTACT US

Creative Industries Centre,  
Glaisher Drive Wolverhampton Science Park,  
Wolverhampton United Kingdom,  
WV10 9TG

**Email queries:**

[info@volmed.co.uk](mailto:info@volmed.co.uk) / [sales@volmed.co.uk](mailto:sales@volmed.co.uk)

**[www.volmed.co.uk](http://www.volmed.co.uk)**