



Serving Healthcare Since 1971



GENERAL SURGICAL INSTRUMENTS CATALOG

www.volmed.co.uk

About Us

With over 50 years of industry experience, VOLMED is a family-owned company specializing in the supply of high-quality surgical instruments to the medical and surgical markets across the United Kingdom, Europe, and globally. Our extensive expertise enables us to assist with project research and development, ensuring that we meet the evolving needs of our clients.

Our philosophy is simple: we treat our customers as we would like to be treated, which places customer satisfaction at the forefront of our priorities. This approach naturally aligns with our commitment to delivering products of the highest quality. Ultimately, our goal is to empower you to provide the best possible care and achieve improved health outcomes by streamlining the process of stocking your practice.

We take pride in the breadth of our product offerings, understanding that sourcing medical tools and surgical instruments from multiple suppliers can be both time-consuming and challenging. To address this, we strive to be your one-stop shop for all your requirements.

Our History

Our story began over 50 years ago within a family deeply rooted in the medical instrument manufacturing industry. The foundation of our thriving business was built by a team of highly ambitious professionals with expertise in producing plastic surgery, dental and surgical instruments.

Inspired by the potential for excellence, our founder established VOLMED, a company dedicated to providing outstanding customer service and high-quality products. Since then, our team of dedicated professionals has worked tirelessly to ensure that customers worldwide receive only the finest surgical instruments and related services. We never compromise on quality, continually investing in our team, resources, equipment, and processes. Every detail is meticulously examined to ensure we can adapt to changing demands while maintaining the highest standards of quality.

Our Promise

At VOLMED, our commitment to excellence is unwavering. We strive to achieve the highest standards in everything we do, working relentlessly to provide exceptional service to our customers and everyone associated with our company. We view it as our responsibility to contribute to a better future for the entire community, and we are dedicated to making this vision a reality.

Quality

Our commitment to quality is reflected in our CE and UKCA-marked instruments, manufactured in accordance with stringent global quality assurance standards, including ISO 9001:2015, ISO 13485:2016, and the Medical Device Directive 93/42/EEC.


Our Mission

We are committed to manufacturing and supplying products that meet the highest standards of strength, accuracy, and reliability. Recognizing the critical role our equipment plays in saving lives and improving health, we take this responsibility seriously. Our mission is to simplify the procurement process by offering a comprehensive range of high-quality surgical products under one roof.

Products

At VOLMED, surgical instruments are at the core of our offerings. We provide a wide range of high-quality surgical products suitable for use in operating rooms, wards, veterinary practices, dentistry, and chiropody. Our instruments are crafted to strict standards using a blend of modern technology and traditional craftsmanship, in full compliance with the Medical Devices Directive 93/42/EEC. Complete traceability is maintained throughout the manufacturing process.

To ensure exceptional product performance and reliability, all our products undergo rigorous engineering analysis and stringent inspection methods. VOLMED instruments are guaranteed for five years against manufacturing defects and are made from premium surgical-grade stainless steel, in accordance with ISO 7153-1:2016 specifications. Explore our product listings and catalogues to discover our comprehensive range of surgical instruments.



Intestinal & Stomach


VOLMED

VOLMED

Serving Healthcare Since 1971


Intestinal & Stomach


Serving Healthcare Since 1971





35

35







Lockwood
V-035-085
20 cm / 8"




Williams
V-035-086
16 cm / 6 1/4"




Littlewood
V-035-190
19 cm / 7 1/2"





Denis-Browne
V-035-300
18 cm / 7"






Young
V-035-290
21 cm / 8 1/4"



Child
V-035-280
24 cm / 9 1/2"



Babcock
Straight



Babcock
Angled 45°


V-035-810	14 cm / 5 1/2"	V-035-800
V-035-281	16 cm / 6 1/4"	V-035-801
V-035-282	18 cm / 7"	V-035-802
V-035-284	20 cm / 8"	V-035-803
V-035-285	24 cm / 9 1/2"	V-035-804
V-035-286	30 cm / 12"	V-035-805
Straight		Angled 45°

Intestinal & Stomach

Tissue and Organ Holding Forceps

Gewebe-und Organ-FaBzangen

Pinzas para tejidos Y organos




13 mm jaws

V-035-287

18 cm / 7"

V-035-289

20 cm / 8"




25 mm jaws

V-035-288


20 cm / 8"

V-035-291

23 cm / 9"



Duval




Duval

V-035-1000

28 cm / 11"

25 mm wide jaws



Allis-Atrauma



V-035-294

16 cm / 6 1/4"



V-035-295

20 cm / 8"




V-035-296

25 cm / 10"



V-035-297

30 cm / 12"



Babcock-Atrauma



V-035-292

16 cm / 6 1/4"



V-035-293


20 cm / 8"



V-035-298

18 cm / 7"


Narrow




V-035-299

20 cm / 8"

Medium




Scheibe



V-035-1010

23 cm / 9"


Wide



V-035-720	15 cm / 6"
V-035-721	18 cm / 7"
V-035-722	20 cm / 8"
V-035-723	23 cm / 9"
V-035-724	25 cm / 10"
V-035-725	30 cm / 12"
3:4 Teeth	

Allis-Baby


V-035-301
13 cm / 5"



V-035-302	15 cm / 6"
V-035-730	18 cm / 7"
V-035-731	20 cm / 8"
V-035-732	23 cm / 9"
V-035-733	25 cm / 10"
V-035-734	30 cm / 12"
4:5 Teeth	


V-035-740	15 cm / 6"
V-035-741	18 cm / 7"
V-035-742	19 cm / 7½"
V-035-743	23 cm / 9"
V-035-744	25 cm / 10"
V-035-745	30 cm / 12"
5:6 Teeth	

Allis




Boys-Allis

V-035-306
15 cm / 6"
7:8 Teeth




Allis-Adair

V-035-307
15 cm / 6"
10:11 Teeth




Judd-Allis

V-035-308 15 cm / 6"
V-035-309 19 cm / 7½"




Lockwood-Allis

V-035-070 16 cm / 6¼"
V-035-080 20 cm / 8"



Allis

V-035-081
20 cm / 8"
5:6 Teeth



Thoms-Allis

V-035-082
20 cm / 8"
6:7 Teeth


Intestinal & Stomach

Tissue and Organ Holding Forceps, Intestinal Clamp Forceps

Gewebe- und Organ-Faßzangen, Darmkelmmen


Pinzas para tejidos Y organos, Pinzas-clamps intestinales

VOLMED




TC Gold
VT-035-351 16 cm / 6¼"
VT-035-380 18 cm / 7"
VT-035-352 20 cm / 8"
VT-035-381 24 cm / 9½"
VT-035-382 30 cm / 12"
 Straight
 Normal profile

Babcock




TC Gold
VT-035-370 16 cm / 6¼"
VT-035-372 18 cm / 7"
VT-035-371 20 cm / 8"
VT-035-373 24 cm / 9½"
VT-035-374 30 cm / 12"
 Angled 45°
 Normal profile

Babcock




Duval
TC Gold
VT-035-353
 20 cm / 8"


Normal profile




Allis
TC Gold
VT-035-400 15 cm / 6"
VT-035-401 19 cm / 7½"
VT-035-402 23 cm / 9"
VT-035-403 25 cm / 10"
VT-035-404 30 cm / 12"




Dannes
V-035-368
 23 cm / 9"




Allen
V-035-1020 20.5 cm / 8¼"
V-035-1025 24 cm / 9½"




Baby-Allen
V-035-371
 16.5 cm / 6½"




V-035-355



V-035-356



Baby-Doyen
17 cm / 6 3/4"



Doyen


V-035-357 21 cm / 8 1/4"
V-035-358 23 cm / 9"




V-035-359
V-035-361




V-035-362 17 cm / 6 3/4"
V-035-363 21 cm / 8 1/4"
V-035-364 24 cm / 9 1/2"




V-035-1030 17 cm / 6 3/4"
V-035-1032 21 cm / 8 1/4"
V-035-1034 24 cm / 9 1/2"




Mayo-Robson / Doyen



V-035-367



V-035-366










Mayo Robson
25 cm / 10"

Intestinal & Stomach

Atrauma Intestinal and Stomach Clamp Forceps with Serrations

Atrauma Darm-und Magnklemmen mit Zahnung

Atrauma Pinzas-Clamps intestinales con dentado






Intestinal & Stomach


Atrauma Intestinal Clamps, Sigmoid Anastomosis Forceps

Atrauma Darmklemmen, Rektum-Anastomosenklemmen

Atrauma Pinzas intestinales, Pinzas para anastomosis rectales




V-035-403




Scudder-Atrauma
33 cm / 13"




V-035-404




Brunner
V-035-405
24 cm / 9½"







Lane
V-035-630
30 cm / 12"



V-035-631



Lane
V-035-640
30 cm / 12"



V-035-641


Intestinal & Stomach

Intestinal and Pylorus Clamps

Darmquetschklemmen

Pinzas intestinales para piloro

VOLMED



Payr-Baby

V-035-650

15 cm / 6"

With pin

Payr

V-035-660

With pin

V-035-661

Without pin

21 cm / 8 1/4"

Payr

V-035-670

With pin

V-035-671

Without pin

29 cm / 11 1/2"

Payr

V-035-680

With pin


V-035-681

Without pin

V-035-679

With groove

035 cm / 14"




DeMartel-Wolfson
Anastomosis Clamps
Sets consist of
2 clamps with short grip
1 clamp with long grip


V-035-682 Set of 3, long clamps
V-035-683 5 cm / 2", clamp with short grip
V-035-684 5 cm / 2", clamp with long grip

V-035-685 Set of 3, long clamps
V-035-686 7.5 cm / 3", clamp with short grip
V-035-687 7.5 cm / 3", clamp with long grip

V-035-688 Set of 3, long clamps
V-035-689 10 cm / 4", clamp with short grip
V-035-691 10 cm / 4", clamp with long grip



V-035-694
26 cm / 10 1/4"
Cosing Forceps



V-035-692 8.5 cm / 3 1/4"
V-035-693 23 cm / 9"
Clamp holder



Serving Healthcare Since 1971

CONTACT US

Creative Industries Centre,
Glaisher Drive Wolverhampton Science Park,
Wolverhampton United Kingdom,
WV10 9TG

Email queries:

info@volmed.co.uk / sales@volmed.co.uk

www.volmed.co.uk